


**REQUIRED SETTING FOR ACCESSING [www.4Medica.com](http://www.4Medica.com)**

---


1. Internet Explorer 5.5 SP2 or 6.0 is required. Updates can be found at [www.microsoft.com](http://www.microsoft.com). Look for the “downloads” link and choose the correct version.
2. Verify that the client does not have any “Pop up Blocker” software installed. Any “Pop up Blocker” software will need to be turned off before each login to 4Medica.
3. Open Internet Explorer and click on “Tools”. Then, click on “Internet Options”.




Click on “Delete Files” and then “OK”. This may take 30 seconds or more to complete.

**REQUIRED SETTING FOR ACCESSING [www.4Medica.com](http://www.4Medica.com)**

4. Contained within the same 'General' tab, click on the 'Settings' button




The 'Check for newer versions of stored pages' option should be set to 'Automatically'




**REQUIRED SETTING FOR ACCESSING [www.4Medica.com](http://www.4Medica.com)**

5. Also under the “Advanced” tab, scroll down to the “Printing” section and verify that there is a check mark next to “Print background colors and images”.




**REQUIRED SETTING FOR ACCESSING [www.4Medica.com](http://www.4Medica.com)**

6. Clients using Windows XP and Internet Explorer 6.0 Service Pack 2 (SP2) will be required to add 4Medica to their trusted sites. Click on the “Security” tab and verify the “Trusted Sites”.



Click on the “Sites” button and add <https://www.4medica.net> or <https://www.4medica.com>. Remove the check mark for “require server verification (https for all sites) in this zone”.





**REQUIRED SETTING FOR ACCESSING [www.4Medica.com](http://www.4Medica.com)**

7. Close the "Internet Options" box and click on "File" and "Page Setup". Remove the "Footer" information and change all the "Margins" settings to "0.5".

