

STOOL COLLECTION CONTAINERS

If request is for:

- Stool Culture
- Enteric Pathogens
- Salmonella/Shigella
- C & S
- Enterohemorrhagic E. coli (e.g., E coli O157)



Use.....

C & S Media. Orange cap.
Store at room temperature

Specimen Container Code **PPCS**

If request is for:

- Stool Parasites
- O & P
- Giardia
- Amoeba
- Intestinal Parasites
- Worms
- Giardia Antigen
- Cryptosporidium Antigen



Use....

Para-Pak. Pink cap. Store at
room temperature

Specimen Container Code

If request is for:

- C. difficile toxin/culture.
Freeze if holding > 24 hrs
- Fecal Leukocytes
- Fecal WBC
- Rotavirus



Use....

Clean Transport Container.
Store at 2-8⁰ (refrigerate)
unless otherwise noted.

Specimen Container Code **ST**



Consolidated Laboratory Services

"Continued Commitment to Quality and Service"

308-5600